| Name of Document | M/S Gulf Marketing Service Rawalpindi | M/S Martin Dow Pvt. Ltd Lahore | M/S A.M Sale & Service Multan | M/S Briogen Pvt. Ltd Lahore |
|---|--|-----------------------------------|----------------------------------|--------------------------------|
| Tender Purchase Receipt Attached | Yes | Yes | Yes | Yes |
| NIC attested copy attached | Yes | Yes | Yes | Yes |
| NTN attested copy attached | Yes | Yes | Yes | Yes |
| STN attested copy attached | Yes | Yes | Yes | Yes |
| One Year Business History | Not Relevant Not Signed | Yes | Yes | Yes |
| Authority Letter / distribution certificate of manufacture / Importer | Yes | Yes | Yes | Yes |
| Establishment License/ License to import | Yes | Yes | Yes | Yes |
| Name of Manufacturer of instrument | ORPHEE | Sweden | Nihon Kohden | France |
| Model Name and Model No. | Mythic TM 18 | Swelab Alfa Plus | Celltac MEK-6510/K | HEMIX 3-60 |
| Country of origin | Switzerland | Sweden | Japan | France |
| Country of manufacturing of instrument | Switzerland | Sweden | Japan | France |
| Quality Certificate | Yes | Yes | Yes | Yes |
| Affidavit regarding non black listing of the firm | Yes | Yes | Yes | Yes |
| Affidavit regarding Price Reasonability Certificate of the firm | Yes | Yes | Yes | Yes |
| TECHNICALLY ELIGIBILITY OF THE FIRM | Not Eligible | Eligible | Eligible | Eligible |
| | | | | |

| | M/S Gulf Marketing Service Rawalpindi | M/S Martin Dow Pvt. Ltd Lahore | M/S A.M Sale & Service Multan | M/S Briogen Pvt. Ltd Lahore |
|---|--|--------------------------------|----------------------------------|--------------------------------|
| Detailed Requirements: | | | | |
| Automated hematology counter | Yes | Yes | Yes | Yes |
| Auto probe cleaning/wiping | Yes | Yes | Yes | Yes |
| RS232/USB interface | Yes | Yes | Yes | Yes |
| 220V, 50 Hz operated | Yes | Yes | Yes | Yes |
| User Adjustable Settings: | | | | |
| Through put 60 samples/hour or more | Yes | Yes | Yes | NO. |
| 3- Parts differential with minimum 20 parameters | Yes | Yes | Yes | Yes |
| Histograms for WBC RBC and platelet | Yes | Yes | Yes | NO. |
| Sample volume Max 120ul of whole | Yes | Yes | Yes | NO. |
| blood or less | | | | |
| Displayed Parameters: | | | | |
| LED displayed with touch screen | Yes | Yes | Yes | Yes |
| Auto calibration programme | Yes | Yes | Yes | Yes |
| Patient data entry and specimen recognition by keyboard/LED touch | Yes | Yes | Yes | Yes |

| screen/Bar Code | Yes | Yes | Yes | Yes |
|--|-----|-----|-----|-----|
| Patient data and results storage of 5,000 results or more | Yes | Yes | Yes | Yes |
| Accessories: | | | | |
| External laser printer | | | | |
| Consumables, reagents, calibrations, controls and start up (free of cost) | | | | |
| Compatible imported online sine wave UPS with battery backup for 30 Minutes (Emerson, Liebert, Chloride, MGE, APC or equivalent) | Yes | Yes | Yes | Yes |
| Operation Manual with a soft copy | | | | |
| Service manual with a soft copy | | | | |
| One Backup model with same specification | | | | |
| AC for temperature maintenance | | | | |
| Free of cost parts replacement, services and maintenance & visits of engineer | | | | |
| Country of origin | | | | |
| (Europe, Japan, USA) | | | | |
| Complete Blood Count (CBC) | Yes | Yes | Yes | Yes |

| DECISION OF TECHNICAL | Non Responsive | Responsive | Responsive | Non Responsive |
|------------------------------|-------------------------------|------------|------------|---------------------------|
| ADVISORY COMMITTEE | | | | |
| | | | | |
| Reason of Rejection | Platelet count and TLC beyond | | | Through put less than |
| | the total allowable error | | | required |
| | the total allowable error | | | Platelet Count and MCV is |
| | | | | beyond total allowable |
| | | | | error |
| | | | | Cleaning time is 10-15 |
| | | | | minutes (too long) after |
| | | | | detection of clot |
| | | | | Sample volume more than |
| | | | | specification |